

New England Biolabs Product Specification

| | |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>Product Name:</i> | <i>PI-PspI</i> |
| <i>Catalog #:</i> | <i>R0695S/L</i> |
| <i>Concentration:</i> | <i>5,000 units/ml</i> |
| <i>Unit Definition:</i> | <i>One unit is defined as the amount of enzyme required to cleave 1 µg of pAKR7 XmnI-linearized Control Plasmid in 1 hour at 65°C in a total reaction volume of 50 µl.</i> |
| <i>Shelf Life:</i> | <i>24 months</i> |
| <i>Storage Temp:</i> | <i>-20°C</i> |
| <i>Storage Conditions:</i> | <i>300 mM NaCl, 10 mM Tris-HCl (pH 7.4), 1 mM DTT, 0.1 mM EDTA, 50% Glycerol, 500 µg/ml BSA</i> |
| <i>Specification Version:</i> | <i>PS-R0695S/L v1.0</i> |
| <i>Effective Date:</i> | <i>10/02/2013</i> |

Assay Name/Specification (minimum release criteria)

Endonuclease Activity (Nicking) - A 50 µl reaction in NEBuffer PI-PspI containing 1 µg of supercoiled PhiX174 DNA and a minimum of 15 Units of PI-PspI incubated for 4 hours at 65°C results in <20% conversion to the nicked form as determined by agarose gel electrophoresis.

Exonuclease Activity (Radioactivity Release) - A 50 µl reaction in NEBuffer PI-PspI containing 1 µg of a mixture of single and double-stranded [³H] *E. coli* DNA and a minimum of 50 units of PI-PspI incubated for 4 hours at 65°C releases <0.1% of the total radioactivity.

Ligation and Recutting (Terminal Integrity) - After a 5-fold over-digestion of pAKR7-XmnI DNA with PI-PspI, ~75% of the DNA fragments can be ligated with T4 DNA ligase in 16 hours at 16°C. Of these ligated fragments, ~75% can be recut with PI-PspI.

Non-Specific DNase Activity (16 Hour) - A 50 µl reaction in NEBuffer PI-PspI containing 1 µg of pAKR7-XmnI DNA and a minimum of 5 Units of PI-PspI incubated for 16 hours at 65°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 10/02/2013

Derek Robinson
Quality Approver

